



M-MLV III Green One-Step qRT-PCR(Low ROX) Kit

产品信息:

组成	MT612-01 100×20μl rxns	MT612-02 500×20μl rxns
2x Green One-Step mix(Low ROX)	1×1ml	5×1ml
M-MLV III Enzyme mix	1×60μl	1×300μl
DEPC-H ₂ O	1×1.8ml	2×1.8ml

Storage and stability:

is shipped on Dry Ice and can be stored for up to 12 months at -20°C, or up to 2 weeks at 4°C. Repeated freeze/thaw cycles should be avoided.

Notes: Research only

Description:

The M-MLV III Green One-Step qRT-PCR(Low ROX) Kit has been formulated for highly reproducible first-strand cDNA synthesis and subsequent real-time PCR in a single tube. A combination of the latest advances in buffer chemistry together with a reverse transcriptase and hot-start DNA polymerase system, ensures that The M-MLV III Green One-Step qRT-PCR(Hi ROX) Kit produces fast, highly-specific and ultra-sensitive one-step RT-qPCR.

The M-MLV III Green One-Step qRT-PCR(Low ROX) Kit consists of a 2x M-MLV III Green One-Step mix, as well as separate reverse transcriptase and RNase Inhibitor.

PCR Reaction Conditions(for a 20μl reaction)

2x Green One-Step mix(Low ROX)	10μl
10μM Forward Primer	0.8μl
10μM Reverse Primer	0.8μl
M-MLV III Enzyme mix	0.6μl
Template and Primers	as required
ddH ₂ O	up to 20μl

PCR cycling conditions:

The following RT-qPCR conditions are suitable for the M-MLV III Green One-Step qRT-PCR(Low ROX) Kit with the majority of amplicons and real-time PCR instruments. However, the cycling conditions can be varied to suit different machine-specific protocols. M-MLV III Green One-Step qRT-PCR(Low ROX) Kit is compatible with either 3-step or 2-step cycling:

3-step cycling:

Step	Temp	Time	Cycles
Reverse transcription	45°C	10min	1
Polymerase activation	95°C	2min	1
Denaturation	95°C	5s	40
Annealing	60°C	10s	
Extension	72°C	5-20s	

2-step cycling:

Step	Temp	Time	Cycles
Reverse transcription	45°C	10min	1
Polymerase activation	95°C	2min	1
Denaturation	95°C	5s	40
Annealing/ Extension	60-65°C	20-30s	

Optional analysis: After the reaction has reached completion, refer to the instrument instructions for the option of melt-profile analysis.

The conditions above are intended as a guide only; conditions will vary from reaction to reaction and may need optimization.

Important considerations and optimization

Primers: The sequence and concentration of the primers, as well as amplicon length, can be critical for specific amplification, yield and overall efficiency of any RT-qPCR. We strongly recommend taking the following points into consideration when designing and running your RT-qPCR:

- use primer-design software, such as Primer3 (<http://frodo.wi.mit.edu/primer3/>) or visual OMP™ (<http://dnasoftware.com/>). Primers should have a melting temperature (T_m) of approximately 60°C
- optimal amplicon length should be 80-200bp, and should not exceed 400bp
- final primer concentration of 400nM is suitable for most SYBR-Green based reactions, however to determine the optimal concentration we recommend titrating in the range 0.1-1µM
- use an equimolar primer concentration
- when possible, use of intron spanning primers to avoid amplification from genomic DNA

Template: It is important that the RNA template is intact and devoid of DNA or contaminating inhibitors of both reverse transcription and PCR. The recommended amount of template for one-step RT-qPCR is dependent upon the type of RNA used.

- **total RNA:** purified total RNA can be used in the range from 1pg to 1µg per 20µl reaction
- **mRNA:** purified mRNA can be used from 0.01pg per 20µl reaction

MgCl₂: The MgCl₂ concentration in the 1x reaction mix is 3mM. In the majority of real-time PCR conditions this is optimal for both the reverse transcriptase and the hot-start DNA polymerase. If necessary, we suggest titrating the MgCl₂ to a maximum of 5mM.

RT-PCR controls: It is important to detect the presence of contaminating DNA that may affect the reliability of the data. Always include a no-RT control, by omitting the reverse transcriptase from the reaction.

Optional ROX: The M-MLV III Green One-Step qRT-PCR Kit is premixed with ROX (5-carboxy-X-rhodamine, succinimidyl ester), so that where necessary, ROX fluorescence can be optionally detected on certain real-time instruments. If your real-time instrument has the capability of using ROX and you wish to use this option, then this option must be selected by the user in the software.

Instrument compatibility

The M-MLV III Green One-Step qRT-PCR Kit has been optimized for use in SYBR Green-based real-time PCR on the real-time instruments listed in the following compatibility table, each of these instruments having the capacity to analyze the real-time PCR data with the passive reference signal either on or off. The kit is also compatible with several instruments that do not require the use of ROX, such as the Qiagen (Corbett) Rotor-Gene™ 6000, the Bio-Rad CFX96 or the Roche LightCycler® 480.

Manufacturer	Model
ABI	7000, 7300, 7700, 7900, 7900HT, 7900HT FAST, StepOne™, StepOne™ Plus

BM180719